



CLEAN BY 4D[®]

Daily cleaner disinfectant



Clean by 4D is a healthcare disinfectant that delivers effective performance and provides a restorative clean. Use Clean by 4D daily as part of a simple and effective cleaning and disinfection program for your entire facility. With a fantastic fresh scent, Clean by 4D leaves the room clean and disinfected!

EPA Registration Number: 5741-30

CLEAN BY 4D®

FEATURES

- Effective against Antibiotic-Resistant Bacteria**
- Kills *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Salmonella enterica*, *Escherichia coli*, *Klebsiella pneumoniae* (NDM-1)
- Sanitizes in 90 seconds
- Kills Hepatitis C Virus (HCV), Hepatitis B Virus (HBV), HIV-1 (AIDS Virus), and Influenza A Virus on hard, nonporous environmental surfaces

** Kills *Staphylococcus aureus*-MRSA, *Enterococcus faecalis*-VRE and *Acinetobacter baumannii* - Ampicillin Resistant

BENEFITS

- One product for all of your cleaning and disinfection needs
- Works great to remove stains and odors from restroom tile and grout
- Quickly breaks down urine to eliminate odors and remove stains
- Provides a restorative clean

RECOMMENDED USE


- Schools, public buildings, restaurants, lodging/hospitality establishments, athletic/recreation facilities, and restrooms
- Toilet bowls and urinals, tub and sink surfaces and fixtures, showers, floors, and other hard, non-porous environmental surfaces

SPECIFICATION DATA

Dilution	2 oz/gal
pH	4.0 – 6.0
Color	Indigo
Fragrance	Light, Fresh
Specific Gravity @ 24°C/75°F	1.031
Stability: Shelf @ 24°C/75°F	1 year minimum
Stability: Freeze/Thaw	Will withstand 3 cycles
EPA Registration Number	5741-30

ADDITIONAL RESOURCES

Scan the QR code to access this product's safety information and documentation.



EFFICACY DATA

Note: This product contains hydrogen peroxide. Do not use this product on brass. Do not use with metal drum pumps, check valves, or equipment.

BACTERICIDAL ACTIVITY

This product is bactericidal on hard inanimate surfaces modified in the presence of 5% organic serum and 200 ppm hard water.

When used as directed for disinfection, this product is effective against *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Salmonella enterica*, *Escherichia coli*, *Klebsiella pneumoniae* (NDM-1), Methicillin Resistant *Staphylococcus aureus*-MRSA, Vancomycin Resistant *Enterococcus faecalis*-VRE, and *Acinetobacter baumannii*.

Leave surfaces wet for 10 minutes.

*VIRUCIDAL ACTIVITY

This product is virucidal in the presence of 5% organic serum and 200 ppm hard water.

When used as directed, this product is effective against Hepatitis C Virus (HCV), Hepatitis B Virus (HBV), HIV-1 (AIDS Virus), and Influenza A Virus.

* Kills HCV, HBV and HIV on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Hepatitis C Virus (HCV), Hepatitis B Virus (HBV) and Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS).

Leave surfaces wet for 10 minutes.

MILDEWSTAT ACTIVITY

When used as directed for sanitization this product is effective against *Aspergillus niger*.


Allow to air dry.

NON-FOOD CONTACT SANITIZING HARD NON-POROUS SURFACES

When used as directed for sanitization this product is effective against *Staphylococcus aureus* and *Klebsiella pneumoniae*. Prepare fresh solution daily or more often if solution becomes visibly dirty, diluted, or soiled.

Treated surface must remain wet for 90 seconds.

ORDERING INFORMATION



SAFETY: For institutional and industrial use only. Be sure to read all directions, precautionary and first aid statements on product labels before using this or any other Spartan product. If questions remain, consult your employer or a physician. Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor and at www.spartanchemical.com. Workplace labels are available in English, Spanish, and French.

GUARANTEE: Spartan's modern manufacturing and laboratory control ensures uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within 12 months of the date of manufacture.